

Comprehensive Mass Spectral Calibration to Achieve High Mass Accuracy and Parameter Free Peak Detection

Yongdong Wang and Ming Gu, Cerno Bioscience, 5 Science Park, Suite 13, New Haven, CT 06511



Introduction


In qualitative MS analysis, there is a growing need for measuring masses at higher and higher mass accuracy (5ppm or even 100ppb) to allow for unambiguous molecular identification in metabolism, biomarker discovery, and proteomics research. In quantitative MS analysis, constant efforts are being made to lower the MS detection limit and increase the quantitative accuracy. Careful examination of MS systems and the overall MS analysis process, however, reveals that the mass accuracy in a given MS system is largely limited by the rudimentary mass calibration schemes used and the detection limit is mostly imposed by the highly heuristic peak detection methods available. This paper will introduce a Highly Accurate Mass Spectral Calibration Approach (HAMSCA) implemented in MSIntegrity™ software to achieve these key objectives in a single-step process (Ref1):

- > High mass accuracy calibration
- > Mass Spectral peak shape correction
- > Noise-filtering without peak distortion
- > Parameter-free and robust peak picking


Methods

Based on the operating principles of various types of MS systems, a comprehensive instrument calibration is performed with an instrument calibration standard giving multiple ions in the mass spectral range of interest. A continuum mass spectral profile data is collected where each known ion represents a corresponding mass shift and peak shape deviation. A calibration function can then be derived to compensate for both the mass shift and the peak shape deviation and applied to successive MS data to convert raw continuum data into fully calibrated continuum data. The fully calibrated MS data can now be subjected to peak detection through elegant mathematical operations without assuming the usual set of peak-picking parameters.


MSIntegrity™ Calibration




Conventional Calibration on Finnigan TSQ Ultra




HAMSCA on Finnigan TSQ Ultra



HAMSCA Signal-to-Noise Gain




Mass Error from Waters qTOF II (Same 587.3153Da as Internal Standard)




Results and Discussion

Peak Centroiding Results from MassLynx™



Peak Centroiding Results from MSIntegrity™




Single and 20-scan Statistics: MSIntegrity™ (AB/Sciex qStar)

Scan Number	158.9646	430.9142	566.8889	702.8638	838.8386	974.8134
340	158.9647	430.9150	566.8880	702.8632	838.8384	974.8140
350	158.9647	430.9148	566.8903	702.8656	838.8366	974.8099
360	158.9649	430.9152	566.8934	702.8629	838.8376	974.8177
370	158.9645	430.9155	566.8862	702.8672	838.8406	974.8143
380	158.9647	430.9139	566.8889	702.8647	838.8384	974.8120
390	158.9641	430.9145	566.8965	702.8673	838.8376	974.8084
400	158.9647	430.9139	566.8889	702.8654	838.8388	974.8101
410	158.9649	430.9140	566.8914	702.8654	838.8384	974.8084
420	158.9648	430.9160	566.8904	702.8670	838.8409	974.8145
430	158.9651	430.9161	566.8918	702.8668	838.8416	974.8159
440	158.9650	430.9153	566.8918	702.8637	838.8387	974.8105
450	158.9647	430.9149	566.8910	702.8669	838.8380	974.8157
460	158.9650	430.9140	566.8912	702.8635	838.8443	974.8150
470	158.9648	430.9151	566.8884	702.8678	838.8407	974.8138
480	158.9648	430.9151	566.8884	702.8678	838.8407	974.8138
490	158.9646	430.9139	566.8886	702.8630	838.8423	974.8140
500	158.9648	430.9147	566.8896	702.8654	838.8384	974.8170
510	158.9643	430.9144	566.8900	702.8630	838.8424	974.8192
520	158.9648	430.9145	566.8896	702.8635	838.8434	974.8099
530	158.9647	430.9139	566.8922	702.8625	838.8387	974.8107
Mean	158.9647	430.9147	566.8889	702.8650	838.8392	974.8131
Single Scan Std Dev (mu)	0.0002	0.0006	0.0015	0.0025	0.0028	0.0032
Single Scan Bias (ppm)	1.5	1.8	1.1	3.7	3.1	3.3
20-scan Bias (ppm)	0.0001	0.0006	0.0009	0.0012	0.0006	-0.0003
20-scan Bias (ppm)	0.6	1.1	1.6	1.7	0.8	-0.3

Single and 20-scan Statistics: Analyst QS™

Scan Number	158.9646	430.9142	566.8889	702.8638	838.8386	974.8134
340	158.9646	430.9139	566.8886	702.8632	838.8443	974.8168
350	158.9642	430.9156	566.8914	702.8639	838.8311	974.8085
360	158.9643	430.9163	566.8943	702.8649	838.8341	974.8212
370	158.9640	430.9161	566.8889	702.8650	838.8348	974.8173
380	158.9642	430.9151	566.8876	702.8651	838.8337	974.8084
390	158.9640	430.9140	566.8829	702.8682	838.8391	974.8173
400	158.9639	430.9151	566.8868	702.8671	838.8369	974.8042
410	158.9639	430.9172	566.8927	702.8692	838.84	974.8206
420	158.9644	430.9168	566.8914	702.874	838.8366	974.8194
430	158.9641	430.9153	566.8946	702.8647	838.8366	974.8164
440	158.9644	430.9145	566.8915	702.8696	838.8386	974.8153
450	158.9644	430.9155	566.8936	702.8691	838.8461	974.8284
460	158.9642	430.9157	566.8905	702.8692	838.8424	974.8053
470	158.9638	430.9141	566.8901	702.8682	838.8382	974.8087
480	158.9639	430.9142	566.8899	702.8713	838.8407	974.8118
490	158.9639	430.9146	566.8895	702.8693	838.8337	974.8152
500	158.9639	430.9147	566.8899	702.8713	838.8407	974.8118
510	158.9641	430.9147	566.8902	702.8682	838.8395	974.8115
520	158.9639	430.9145	566.8895	702.8693	838.8318	974.8102
Mean	158.9640	430.9153	566.8905	702.8680	838.8364	974.8076
Single Scan Std Dev (mu)	0.0002	0.0009	0.0014	0.0017	0.0044	0.0136
Single Scan Bias (ppm)	1.6	2.4	1.1	5.2	5.1	13.8
20-scan Bias (ppm)	-0.0008	0.0011	0.0015	0.0022	-0.0022	-0.0098
20-scan Bias (ppm)	0.6	2.4	2.7	3.2	2.4	-6.8


Loop Injection on API 4000 QTRAP




Mass Accuracy on API 4000 QTRAP (285.1425 and 609.2812Da for Internal Cal)

Injection #	285.1425	388.2055	472.3124	477.2325	508.2659	609.2811	716.4482
1	285.1427	388.2066	472.3124	477.2322	508.2658	609.2811	716.4482
2	285.1428	388.2058	472.3124	477.2326	508.2659	609.2812	716.4480
3	285.1428	388.2078	472.3120	477.2304	508.2658	609.2811	716.4484
4	285.1430	388.2082	472.3125	477.2349	508.2654	609.2811	716.4377
5	285.1428	388.2088	472.3124	477.2326	508.2658	609.2811	716.4486
6	285.1428	388.2078	472.3124	477.2326	508.2658	609.2811	716.4486
7	285.1431	388.2075	472.3094	477.2303	508.2652	609.2814	716.4458
8	285.1428	388.2065	472.3124	477.2327	508.2652	609.2813	716.4454
Mean	285.1434	388.2079	472.3174	477.2306	508.2675	609.2812	716.4411
10-Injection Std Dev (mu)	0.0003	0.0007	0.0004	0.0004	0.0005	0.0001	0.0004


Mass Accuracy Comparison to Analyst™ (AB/Sciex API 4000 QTrap LC/MS)



Raw Data Before Alignment



HAMSCA Mass Spectral Alignment



Conclusions

- > Calibrating for mass alone would not yield the highest possible mass accuracy.
- > Peak shape correction is as important as mass calibration for mass accuracy.
- > With the use of HAMSCA in MSIntegrity, the mass error on a unit mass resolution system can be reduced from 0.1-0.2Da to several mDa with a gain factor approaching 100x.
- > Effective noise filtering of over 3x can be achieved without the usually detrimental peak shape distortion.
- > Unambiguous and parameter-free peak picking preserves all quantitative information such as isotope ratios while taking advantage of all available ion signals.
- > HAMSCA also enables quantitative mass spectral alignment, a key requirement in applying multivariate statistics to mass spectral data.

Notes

1. Ming Gu, Yongdong Wang, and Don Kuehl, *Spectroscopy (Current Trends in Mass Spectrometry)*, 2005, p22-26.