

Formula Determination and Relative Quantitation from Overlapping Ion Signals

Yongdong Wang and Ming Gu
Cerno Bioscience, Danbury, CT, USA

Introduction

Comparing to other analytical technologies, mass spectrometry is typically considered of high selectivity, due to its relatively high resolving power, even at the conventional unit mass resolution. Overlapping mass spectral signals, however, can be observed frequently in mass spectrometry, even on higher resolution MS systems. Some examples include:

- The overlap of the two different forms of the same drug metabolites, one with and one without isotope labeling (^{14}C or ^{15}N), resulting in ions offsetting each other by the exact mass of one mass unit with very small mass defects, causing interferences to the two ions involved.
- The fragment ions from EI GC/MS or LC/MS/MS experiments.
- Product oxidation or reduction through the loss or gain of H_2 . This paper will present a novel approach to both identify and quantify these mutually overlapping ions.


Methods

While it is sometimes possible to move to MS systems of higher and higher resolving power to separate these overlapping mass spectral signals, a different approach will be presented in this paper to mathematically resolve these mutually interfering signals. To accurately deconvolute these overlapping signals, a unique mass spectral calibration has to be performed that calibrates not only m/z axis but, more importantly, mass spectral peak shape. With the Spectral Accuracy achieved through this calibration, the overlapping mass spectral signals can now be resolved reliably and mathematically, leading to both qualitative (e.g., elemental compositions) and quantitative results (e.g., the relative concentration of the overlapping ions), even at unit mass resolution.


Theory

Unlike the conventional MS calibration, a novel calibration involving both m/z and MS peak shape allows for exact isotope modeling and formula determination on even a unit mass resolution system, provided that the mass spectral signal is free from significant interferences arising from coexisting ions.

MassWorks CLIPS Formula ID




In the presence of mass spectral interferences, the Spectral Accuracy and formula ID will be compromised. When the mass spectral interference is from an ion related to the ion of interest, however, the same exact isotope modeling can be expanded to include the interference ions in the Spectral Accuracy and formula ID process, as shown by the example below.






Results and Discussion

EI GC/MS




In GC/MS experiments with EI source, the electron beam can often strike a gas molecule hard enough for it to lose one or more hydrogen atom(s) during the electron impact process. This can happen to both the molecular ion itself or its fragments, especially if the molecular ion is a hydrocarbon. In the above graph, the molecular ion loses both one and two hydrogen atoms, forming a three component mixture in the observed mass spectrum. With the mass spectral calibration performed, all these three components can be mathematically resolved down to the random noise level, allowing for both qualitative formula determination and relative quantitation of all three components involved.

LC/MS/MS Fragments ID




In this LC/MS/MS experiment with a triple quad from Waters Quattro, the MassWorks calibration is performed using known fragments ions coming from the same precursor ion to achieve a calibration on Q3 which serves as the mass analyzer. Alternatively, one could have used the known fragments of an entirely different precursor ion. The various fragments may or may not overlap with each other but could be separated mathematically for structural elucidation purpose. The isolation window for the MS/MS experiment is set to be open, allowing for all isotopes of the precursor ion to come through for this analysis.

HiRes Quad w/Pharmaceuticals - A Self Calibration Case Study




For drug compound reserpine, an impurity is typically observed whose m/z is close enough to become an interference to the standard reference line at 609.2812Da. While it would have required a calibration process in order to perform both identification and relative quantitation of this impurity, we tested a self calibration process called sCLIPS using higher resolution quadrupole MS data available on a Thermo Quantum instrument here with interesting and useful results. Due to the lack of mass calibration, a large mass tolerance of 100mDa is used to come up with a list of 488 formula candidates. With the peak shape calibration, however, the high Spectral Accuracy of 97.8% allows the correct formula to be ranked as the 22nd hit with unbiased relative quantitative information (9.7% impurity). The results point to possible oxidation via the loss of two hydrogen atoms.

Formula	Monoisotopic Mass	Mass Error (mDa)	Relative Intensity (%)	Spectral Accuracy (%)	Rank	
1	C ₂₂ H ₂₇ O ₅	-80.0779	100.0000	98.3322	3.346	22.0
2	C ₂₂ H ₂₇ O ₅	-74.2600	122.2101	98.3113	3.386	14.0
3	C ₂₂ H ₂₇ O ₅	-48.8465	-73.8208	98.2501	3.505	22.0
4	C ₂₂ H ₂₇ O ₅	-4.4500	-13.5538	98.2252	3.537	21.0
5	C ₂₂ H ₂₇ O ₅	2.7744	4.6960	98.2346	3.542	21.0
6	C ₂₂ H ₂₇ O ₅	18.4414	16.7000	98.2354	3.544	20.0
7	C ₂₂ H ₂₇ O ₅	-18.6024	-32.4306	98.1460	3.722	21.0
8	C ₂₂ H ₂₇ O ₅	-14.0138	-23.8662	98.1379	3.726	21.0
9	C ₂₂ H ₂₇ O ₅	-79.8608	-131.5719	98.1148	3.736	20.0
10	C ₂₂ H ₂₇ O ₅	-1478	-1133721	98.0704	3.820	20.0
11	C ₂₂ H ₂₇ O ₅	17.8255	32.2963	98.0699	3.830	16.0
12	C ₂₂ H ₂₇ O ₅	-3.8111	-56.3528	98.0667	3.841	22.0
13	C ₂₂ H ₂₇ O ₅	18.8748	19.2037	98.0663	3.840	14.0
14	C ₂₂ H ₂₇ O ₅	-82.4834	-152.2742	98.0567	3.903	20.0
15	C ₂₂ H ₂₇ O ₅	-61.1203	-102.0771	98.0488	3.915	20.0
16	C ₂₂ H ₂₇ O ₅	53.9833	82.9006	98.0396	3.936	20.0
17	C ₂₂ H ₂₇ O ₅	1.6187	10.4044	97.9954	4.442	20.0
18	C ₂₂ H ₂₇ O ₅	-57.4258	-24.5723	97.8620	4.240	20.0
19	C ₂₂ H ₂₇ O ₅	15.2411	41.6568	97.8684	4.052	21.0
20	C ₂₂ H ₂₇ O ₅	39.5589	64.4608	97.8472	4.119	20.0
21	C ₂₂ H ₂₇ O ₅	72.8503	119.6252	97.8465	4.252	20.0
22	C ₂₂ H ₂₇ O ₅	-68.9569	-80.7881	97.8211	4.267	16.0
23	C ₂₂ H ₂₇ O ₅	-23.2771	-38.8303	97.8088	4.312	20.0
24	C ₂₂ H ₂₇ O ₅	-41.2300	-133.7744	97.8448	4.332	20.0
25	C ₂₂ H ₂₇ O ₅	21.7085	45.9910	97.8362	4.398	20.0
26	C ₂₂ H ₂₇ O ₅	-12.4603	-34.6162	97.8279	4.398	21.0
27	C ₂₂ H ₂₇ O ₅	-46.1055	-76.0725	97.7904	4.415	20.0
28	C ₂₂ H ₂₇ O ₅	-32.2737	-53.1503	97.7791	4.474	21.0
29	C ₂₂ H ₂₇ O ₅	-8.8507	-16.1507	97.7686	4.477	21.0
30	C ₂₂ H ₂₇ O ₅	15.6441	154.6461	97.7685	4.477	20.0



DART AccuTOF Food Analysis



Conclusions

- Though the presence of overlapping ion signals may be difficult or even impossible to overcome through either chromatography or higher resolution MS, there exists an elegant mathematical solution.
- The mathematical solution provides both qualitative results (formulae) as well as quantitative results (relative concentrations).
- At unit mass resolution, this mathematical solution calls for a comprehensive MS calibration involving both m/z and peak shape through known standards.
- At higher-than-unit mass resolution, it is possible to perform both qualitative and quantitative analysis through self calibration, i.e., without the use of any known standards. While a larger mass tolerance window may be needed, which may compromise the performance, the differentiation power of Spectral Accuracy can still produce meaningful results.

Acknowledgments

- Zhe-Ming Gu & Mark Ma (XenoBiotic Lab) for QTrap 4000 LC/MS isotope labeling data.
- James Mullis and Fenghe Qiu (Boehringer Ingelheim, USA) for Agilent GC/MSD EI data.
- Bob Strife, Michele Mangels & Jason Price from P&G (Cincinnati, OH) for the Waters Quattro LC/MS/MS data.
- Leo Okkerse (formerly Thermo Electron) for the high resolution Quantum Max triple quadrupole data of reserpine.
- Tommy Lewander (LSM Lab, Sweden) and Jan Nordin (Chemalys, Sweden) for DART AccuTOF data of olive oils.