

Accurate Mass Measurements of Product Ions for Metabolite Identification on Unit Mass Resolution Mass Spectrometers

Jianyao Wang, Wyeth Pharmaceuticals, 500 Arcola Road, Collegeville, PA
Ming Gu and Yongdong Wang, Cerno Bioscience, Monmouth Junction, NJ



Overview

- Accurate Mass measurements of product ions made available on unit mass resolution mass spectrometers
- Novel algorithms to calibrate mass spectrometers for both mass and peak shape
- Demonstrated data including infusion of drug standards and LC/MS of rat liver microsomal incubations of verapamil

Introduction

Accurate mass (AM) measurements are routinely performed on high resolution mass spectrometers for metabolite identification. However, these instruments are very expensive and often require dedicated and skillful operators. We have developed MSIntegrity™ technologies that allow AM analysis to be conducted on unit mass resolution mass spectrometers¹. Although AM analysis of verapamil metabolites has already been demonstrated previously in full MS mode, further elucidation of the structures of these metabolites requires the recognition of their corresponding fragment ions from collision-induced dissociation experiments. Accurate measurements of these ions are desirable to pinpoint their structures². As an extension to the MSIntegrity technologies, we report here the AM measurements of fragment ions to facilitate metabolite identification on unit mass resolution mass spectrometers.


Methods

The product ion spectra of a set of small pharmaceuticals were acquired from an ABI/Sciex 4000 Q TRAP by infusion, while LC/MS/MS measurements of rat liver microsomal incubation of verapamil were performed on an ABI/Sciex 3200 Q TRAP.

Methods


All mass spectral data were collected in profile mode and processed using the MSIntegrity technologies implemented in the MassWorks™ software. A few or a single known fragment ion(s) from the MS/MS spectra by infusion was selected as calibration standard(s) to internally calibrate the spectra. Some of the ions in the same spectra were tested for mass accuracy. LC/MS/MS fragment ions were generated on-the-fly by information dependent scans. Accurate mass measurements on demethylated verapamil were achieved by applying an external calibration made by the fragments generated from verapamil.

Calibration Procedures




Results and Discussion

The MS/MS Peak Shape before and after Calibration



MS/MS Spectrum of Reserpine




AM of Product Ions of Reserpine (Internal calibration)

Four point calibration					
Name of Ions	Exact Mass	Formula of Ions	Accurate Mass	Delta Mass (mDa)	
Reserpine	606.2812	C ₃₃ H ₄₁ N ₂ O ₉	Call Ion		
Frag1	448.1971	C ₂₂ H ₂₅ N ₂ O ₆	448.1968	-3.4	
Frag2	360.1472	C ₁₈ H ₂₁ N ₂ O ₄	Call Ion		
Frag3	198.0857	C ₁₀ H ₁₁ O ₄	Call Ion		
Frag4	174.0913	C ₁₁ H ₁₃ O ₃	Call Ion		

Two point calibration					
Name of Ions	Exact Mass	Formula of Ions	Accurate Mass	Delta Mass (mDa)	
Frag2	360.1472	C ₁₈ H ₂₁ N ₂ O ₄	Call Ion		
Frag3	198.0857	C ₁₀ H ₁₁ O ₄	198.0830	-2.7	
Frag4	174.0913	C ₁₁ H ₁₃ O ₃	Call Ion		

One point calibration					
Name of Ions	Exact Mass	Formula of Ions	Accurate Mass	Delta Mass (mDa)	
Frag3	198.0857	C ₁₀ H ₁₁ O ₄	Call Ion		
Frag4	174.0913	C ₁₁ H ₁₃ O ₃	174.0924	+1.6	

MS/MS Spectrum of Clozapine




AM of Product Ions of Clozapine (Internal Calibration)

One point calibration					
Name of Ions	Exact Mass	Formula of Ions	Accurate Mass	Delta Mass (mDa)	
Reserpine	606.2812	C ₃₃ H ₄₁ N ₂ O ₉	606.2809	-3.3	
Frag1	268.0963	C ₁₇ H ₁₅ N ₂	Call Ion		
Frag2	270.0769	C ₁₅ H ₁₃ N ₂ O	270.0498	-36	


Note: Large mass errors of m/z of 270 were due to saturation.

Results and Discussion


Saturated and Theoretical Isotope Distribution of C₁₅H₁₃CIN₃




TIC of Verapamil Metabolites




MS/MS of Verapamil



MS/MS of Demethylated Verapamil



Mass Errors of Four Calibration Ions



AM of MSMS of Demethylated Verapamil (External Calibration)

Name of Ions	Exact Mass	Formula of Ions	Accurate Mass	Delta Mass (mDa)
Demethylated Verapamil	441.2753	C ₂₆ H ₃₇ N ₂ O ₄	441.5840	-0.6
Frag1	303.2073	C ₁₈ H ₂₇ N ₂ O ₂	303.4159	5.5
Frag2	289.1916	C ₁₇ H ₂₅ N ₂ O ₂	289.3933	-5.6
Frag3	291.2073	C ₁₇ H ₂₇ N ₂ O ₂	291.4099	-6

Conclusions

- Accurate mass measurements of product ions with 5 ppm can be achieved on unit mass resolution instruments through MassWorks
- Mass spectrometry calibration including both mass and peak shape is critical for high mass accuracy.
- Internal calibration can help achieve highest possible mass accuracy.
- Peak saturation is a key factor causing poor mass accuracy in a quadrupole system.

References

- 1) Gu, M. et al *RCMS*, 2006, 20, 764.
- 2) Zhang, M. et al *JMS*, 2005, 40, 1017.