

Confident Unknown Identification of SVOC Compounds by Combining NIST Library Search with Elemental Composition Determination

Jianping Chen¹, O. David Sparkman², and Ming Gu³

¹Connecticut State Department of Environmental Protection, Windsor, CT; ²University of the Pacific, Stockton, CA; ³Cerno Bioscience, Danbury, CT

Overview

➤ A new method based on mass spectral peak shape calibration was developed for the identification of unknown Semivolatile Organic Compounds (SVOCs).

➤ On a unit mass resolution quadrupole GC/MS instrument, 122 SVOC compounds were identified and/or confirmed by NIST library search and MassWorks formula determination.

➤ Complementary approaches to SVOC unknown ID by NIST Mass Spectra Database and CLIPS search were discussed.

Introduction

Using GC/MS for the rapid identification of releases of hazardous or toxic chemicals requires confirming the analyte identity by a second method. High resolution MS such as TOF may provide this confirmation, but it is not an option due to high cost and limited space available in this mobile lab. Alternatively, a novel method which provides a primary and secondary confirmation with a single analysis using a well established EPA procedure for SVOCs has been developed. The method is based on formula determination by advanced data processing of integer-mass-only results produced by the GC-MS in conjunction with a conventional library search. This dual process dramatically improves the accuracy of chemical identification.

Methods

➤ **Data Acquisition:** All data were acquired in an Agilent 6890-5973N GC/MS system in a raw (profile) scanning mode and followed by NIST MS Database Search and elemental composition determination by MassWorks. A total of 122 SVOC compounds were investigated. About 0.1 mg of each compound were injected and separated by GC/MS according to EPA methods. At the end of run, PFTBA spectra were acquired for 20 seconds to built calibration by MassWorks. The system can hold the calibration for days if not weeks before recalibration.


➤ **NIST MS Database Search:** By an Agilent ChemStation data processing script, the profile spectra were automatically converted into centroid spectra for NIST MS Database Search.

➤ **MassWorks processing:** MassWorks calibrates both the mass position and the mass spectral peak shape function, a key for achieving high mass accuracy. When applied to raw mass spectral data, the raw mass spectrum can be transformed into its calibrated version with mass spectral peaks located at accurate mass positions. Furthermore, the mass spectral peak shape would also be transformed in the same process to a mathematically definable function, a key for achieving high spectral accuracy and CLIPS formula ID.

➤ **Performance Evaluation:** Among the 122 SVOC, 13 compounds with a mass range from 93 to 276 were used as routine performance check before running unknown substance identification. High spectral accuracy of 98.5% and top 3 search ranking are required performance criteria for those 13 compounds before applying it to identify real unknown compounds.


Results and Discussion

MassWorks Calibration+CLIPS Formula ID



For each known ion fragments from the PFTBA standard, both mass and peak shape are adjusted as part of the calibration to enable highly accurate mass determination at unit mass resolution in the presence of isotope interferences, as shown from the $C_6F_5^+$ fragment ion (accurate mass 130.9907 vs. exact mass 130.9920Da).

Fig 1. Mass and Peak Shape Calibration with PFTBA



This calibration is then applied to the data file of GC/MS of SVOC. First, the 13 selected performance checking standards were analyzed. As shown in Table 2, all the standards achieved better than expected spectral accuracy of 98.5% and expected ranking of top 3 on search list. As an example shown in Fig 2, Acenaphthylene ranked as top 1 has spectra accuracy of -0.3 mD and spectral accuracy 99.7%. This great performance is attributed to the novel calibration technology that corrects unsymmetrical peak shape (Fig. 2a) and provides effective isotope pattern matching (Fig 2b) between calibrated and theoretically calculated spectra.

Secondly, accurate mass measurements and

Table 1. CLIPS Search Parameters

Accurate Mass	130.99	Element	Min	Max
Charge	1	C	0	50
Mass Tolerance (mDa)	20	H	0	100
Electron State	CGSE	N	0	15
DBE Min	-10	O	0	16
DBE Max	50	S	0	10
Profile Mass Start (Da)	0.5	F	0	12
Profile Mass End (Da)	4.5	P	0	12
		Br	0	10
		Cl	0	12

Table 2. Mass and Spectral Accuracy for Performance Checking Standards

Target Compounds	Formula	Retention Time	m/z	Rank	Spectral Accuracy	Mass Accuracy	Total Formu
Aniline	C6H7N	4.52	93	1	99.18	-4.1	27
1,3-Dichlorobenzene	C6H4Cl2	5.07	156	1	99.90	2.8	196
Benzyl alcohol	C7H8O	5.45	108	2	99.66	-7.4	62
1,4-Methylphenol	C7H8O	5.51	108	1	99.85	-9.9	62
Ch. Voluolene	C7H8	5.78	107	1	99.91	-5.3	39
N-Nitrosopyrrolidine	C5H10NO	6.18	114	1	99.75	-7.3	72
2,4-Dichlorophenol	C6H4Cl2O	6.36	155	1	99.61	-8.1	292
Naphthalene	C10H8	6.48	128	1	97.45	0.2	116
2,4,5-Trichlorophenol	C6HCl3O	7.94	196	2	98.83	-8.3	479
Acenaphthylene	C12H8	8.46	152	1	99.74	-8.3	229
1,1-Dichlorobenzene	C6H4Cl2	11.39	252	2	98.86	-1.8	273
Benzofluorene	C20H12	11.91	252	1	99.73	-4.4	1469
Benzofluorene	C20H12	16.97	276	1	99.49	-3.7	2192

Fig 2. Raw, Calibrated, and Calculated Spectra of m/z 152





Table 3. Top 10 of Search for m/z 152

Row	Formula	Mono Isotope	Mass Error (mDa)	Mass Error (PPM)	Spectral Accuracy	DBE
1	C12H8	152.0626	-0.3184	-2.0936	99.7377	9
2	C8H6F2	152.0438	18.5263	121.8271	98.9747	6
3	C8H6F2	152.0637	-1.4612	-9.6093	98.9539	6
4	C8H4F4	152.0795	-13.205	-86.8396	98.6873	4
5	C8H6Cl2	152.0473	14.9376	98.2446	98.6661	5
6	C8H6Cl2	152.0726	-12.0949	-79.4829	98.5036	5
7	C7H8N2	152.0503	11.8471	78.5671	98.564	5
8	C7H8N2	152.0548	3.7044	24.311	98.5058	5
9	C8H6F4	152.0449	12.4575	81.9238	98.4455	6
10	C8H4F4	152.068	-3.8892	-24.2609	98.3343	4

formula identification were conducted on all the SVOC. Formula search conditions listed on the Table 1 include a wide range of elements and very open lower and upper limits. The electron state of "odd" electron was selected to focus on the analysis of the compounds having molecular ions which accounts for 91% of the SVOC under investigation. Fig 3. shows the ranking distribution of both NIST MS Database and CLIPS search. 94% of the SVOC appeared on top 3 from CLIPS search lists including 76% 1st, 16% 2nd, 2% 3rd hits for the compounds with molecular ions. NIST library search generated positive ID for 92% of SVOC regardless the

Fig 3. NIST Library and CLIPS Search Rank Distribution



existence of the molecular ions. Poor results from NIST library search most likely were due to complex mass spectra caused by co-eluting compounds. For example, the major peaks of the spectra of N-Nitrosopyrrolidine should be m/z 100 as molecular ion and a dominant fragment at m/z 41 (Fig 4b). However, this compound was co-eluted with other molecules and its spectra appears to be much more complicated (Fig 4a). This resulted in poor match and reverse match values at 392 and 593 respectively by NIST MS Database Search. Instead of matching multiple ions fragmentation pattern by a NIST MS Database Search, MassWorks performs formula determination by analyzing individual ions and is able to deal with co-eluted problems as long as the interesting ions are not overlapped by another unknown ion.

Fig. 5 shows that the elemental composition of N-Nitrosopyrrolidine was found to be the 2nd hit by CLIPS search with mass accuracy of 4.5 mD and spectral accuracy of 99.3%

Fig 4. Complex and pure spectra of N-Nitrosopyrrolidine





Fig 5. Isotope Pattern Match for N-Nitrosopyrrolidine



Row	Formula	Mono Isotope	Mass Error (mDa)	Mass Error (PPM)	Spectral Accuracy	RMSE	DBE
1	C5H8NO	100.0521	-6.0709	-60.6558	99.319	835	3
2	C5H8NO	100.0521	-6.0709	-60.6558	99.319	835	3
3	C5H8NO	100.0688	9.7285	97.2274	98.9866	1,247	2
4	C5H8NO	100.0497	-9.3558	-93.9311	98.4509	1,308	3
5	C5H8NO	100.05	-9.1152	-91.0979	98.4372	1,623	1

Conclusions

➤ With PFTBA calibration performed to each GC/MS run, mass accuracy better than 10mDa and Spectral Accuracy better than 99% can be achieved for reliable formula determination.

➤ Based on both NIST library and CLIPS search, confident identification of unknown SVOCs can be achieved.

➤ When NIST library search results in multiple compounds with similar Match Value (MV) and Reverse Match Value (RMV), MassWorks can provide unique formula determination for unknowns.

➤ NIST MS Database Search can provide positive identification regardless the existence of molecular ions.

References

1. Tyler, A. N. et al. *Anal. Chem. Exact Mass Measurement of Polar Organic Molecules at Low Resolution Using Electrospray Ionization and a Quadrupole Mass Spectrometer*, 1996, 68, 3561.
2. Fiehn, O. et al. *Anal. Chem. Identification of Uncommon Plant Metabolites Based on Calculation of Elemental Compositions Using Gas Chromatography and Quadrupole Mass Spectrometry*, 2000, 72, 3573.
3. Gu, M et al. *RCM Accurate mass filtering of ion chromatograms for metabolite identification using a unit mass resolution liquid chromatography/mass spectrometry system*, 2006, 20, 764.
4. Wang, Y.; Prest, H. *Chromatography (Japan) Accurate Mass Measurement on Real Chromatographic Time Scale with a Single Quadrupole Mass Spectrometer*, 2006, 27, 135.