

Automated Molecular Formula Determination in Open Access LC/MS in Drug Discovery.

Vladimír Čápka¹, Sharon Tentarelli¹, Hongliang (Leo) Xu²

¹AstraZeneca R&D Boston, 35 Gatehouse Dr., Waltham, MA 02451, ²Cerno Bioscience, 14 Commerce Dr., Danbury, CT 06810



Overview

- Molecular formula determination based on mass and spectral isotope pattern can play an important role in compound identification.
- Sample analysis using a single quadrupole MS instrument and MassWorks™ software with isotopic pattern calibration can yield molecular formula determination.
- User-friendly automated software processing has been developed to allow for integration into open-access systems in the drug discovery environment.
- Capabilities and limitations of open-access molecular formula determination are compared.


Introduction

- In drug discovery applications, molecular formula determination based on accurate mass and spectral isotope pattern is typically measured using high resolution and performed by skilled analytical chemists.
- Use of high resolution instrumentation such as TOF in open access is problematic due to required user expertise and operation complexity.
- A typical open-access environment contains single quadrupole MS instruments operated by synthesis chemists with moderate MS expertise.
- MassWorks™ software with isotopic pattern calibration can be used for molecular formula determination on a single quadrupole instrument.
- Interfacing MassWorks™ software with MassLynx™ software can automate much of the calibration and sample analysis process, minimizing the need for end-user expertise.
- Effect of user input of target ion, element list, and target molecular formula was investigated.
- 28 model compounds (MW 292-652 Da) with typical elemental composition in drug discovery were used.

Results

Three Typical Workflows

- Chemist desires confirmation of the presence of a known compound.
- Chemist has limited molecular knowledge about an unknown sample.
- Chemist has a completely unknown sample.
- Chemist inputs expected molecular formula during sample login.
- Chemist inputs target ion of interest and/or elements likely to be present during sample login.
- Report shows confirmation, or failure and proposed alternate formulas.
- Chemist inputs no molecular information during sample login.
- Report shows proposed formulas for predominant ions within spectra based on default element list



Formula Identification Using a Variety of Different Workflows

Table I. Effectiveness of sample identification for 28 test compounds based on different user input during sample login..

Formula	Input → Formula			No input			Target ion			Elements		Elements and target ion				
	Mass Error (mDa)	SA (%)	Rank	Mass Error (mDa)	SA (%)	Rank	Mass Error (mDa)	SA (%)	Rank	Mass Error (mDa)	SA (%)	Rank	Mass Error (mDa)	SA (%)		
C ₂₆ H ₂₈ N ₄ O ₂	8.3	98.8	NF							1	-10.6	98.44	2	0	98.7	
C ₃₀ H ₂₇ ClN ₃ O	-7.4	98.7	TNA										2	-4.5	99.0	
C ₁₅ H ₁₀ N ₄ OS	5.1	99.0	1	-11.3	99.2	1	-11.3	99.2	2	9.3	99.29	3	-9.3	98.6		
C ₁₇ H ₁₈ FN ₄ S	5.2	98.3	TNA										2	0.5	98.7	
C ₁₉ H ₂₃ NO ₂	8.5	99.3	TNA							1	1.7	99.3	TNA	1	5.3	98.2
C ₁₈ H ₂₄ ClN ₃ O ₃	2.7	98.0	3	-14.0	98.1	2	-2.7	98.3	2	-4.0	98.59	1	-6.5	99.0		
C ₁₃ H ₁₆ N ₄ O ₂ S	6.7	98.2	1	-11.5	99.1	NF				3	2.6	98.96	1	2.3	98.7	
C ₁₅ H ₁₇ Cl ₂ NO ₂	9	98.7	1	-7.2	99.1	2	-5.2	99.4	1	7.6	99.22	1	-2.4	99.2		
C ₂₀ H ₂₃ ClN ₃ O	7.1	98.3	NF							3	-15.7	99.0	2	0.8	98.78	
C ₁₉ H ₂₂ ClN ₃ O	15	94.2	TNA												NF	
C ₂₂ H ₂₆ N ₂ O ₂	2.4	99.1	TNA							3	-5.6	99.5	TNA	2	-4.9	98.9
C ₂₃ H ₂₈ N ₂ O ₂	-4.5	99.4	NF							2	-7.6	99.4	2	-3.1	99.11	
C ₂₅ H ₂₄ N ₄ O ₆	3.4	98.5	NF							1	-11.6	97.4	NF	2	-3.9	99.1
C ₂₇ H ₂₈ N ₄ O ₃	-13.7	98.9	1	-3.2	98.9	NF				1	-7.4	99.01	1	-11.1	99.1	
C ₂₆ H ₂₆ N ₄ O ₃	-13.5	98.6	NF							1	-10.2	98.98	1	0.2	99.2	
C ₂₇ H ₂₈ N ₆ O ₃	-12.9	98.7	TNA										2	-1.7	98.5	
C ₂₈ H ₃₁ N ₃ O ₃	-13.9	98.1	TNA										3	-14.1	97.5	
C ₂₇ H ₂₂ F ₂ N ₄ O ₄	-12.6	98.0	TNA							1	-9.5	98.4	TNA	2	-5.4	99.2
C ₂₆ H ₁₉ F ₃ N ₆ O ₂	11.8	97.8	NF							1	11.2	98.81	NF			
C ₂₆ H ₂₇ N ₅ O ₂	-2.9	98.1	NF							2	-5.8	99.17	2	-7.4	98.8	
C ₂₆ H ₂₁ ClN ₆ O ₂ S	-12.1	97.8	NF							1	1.0	98.2	3	10.2	98.39	
C ₂₇ H ₂₃ FN ₆ O ₄	-13.2	98.6	TNA											1	-1.4	99.3
C ₃₀ H ₃₃ N ₅ O ₃	-17	97.9	NF										2	-15.8	95.0	
C ₂₅ H ₂₀ ClN ₅ O ₄	-12.6	99.0	NF							1	-12.5	98.8	NF	2	-8.3	99.0
C ₂₄ H ₁₉ ClN ₄ O ₂	-14.1	98.1	NF							1	-12.4	98.79	1	-2.6	99.0	
C ₂₉ H ₂₃ F ₃ N ₉ O ₂ S	1.5	97.8	NF												NF	
C ₂₈ H ₃₂ ClN ₃ O ₈	-11.4	99.0	TNA											2	1.4	99.0
C ₂₉ H ₃₄ ClN ₃ O ₆	-2.9	99.0	TNA											2	8.0	99.3
Total correctly identified		27			5			11			13				24	

One failed sample (SA=94.2) due to poor ionization.

Unlikely workflow; some sample information is usually known (elements, mass).

Most likely workflow; mass and elements are known from other open access analyses and reaction scheme.

NF = Not Found; correct formula was not within top 3 formulas.
TNA = Target Not Analyzed; target ion was not the largest ion within the spectrum.

Materials and Methods

Equipment Setup	
Instrument	Waters ZQ single quadrupole Agilent 1100 HPLC
Column	Waters XBridge™ C18, 3.5 μm, 2.1 x 20 mm
Flow Rate	1.8 mL/min
Injection Volume	5 μL
Mobile Phase	A: 0.1% formic acid in H ₂ O B: 0.1% formic acid in CH ₃ CN
Gradient	5 to 95% B in 1.75 min
MS detection, ESI+	Continuum mode, 175 – 695 amu at 350 amu/s
Resolution	Unit, ~0.5 amu peak width half height
Calibrant	terbutaline ([M+H] ⁺ 226.1443), 100 μg/mL aspartame ([M+H] ⁺ 295.1294), 1000 μg/mL buspirone ([M+H] ⁺ 386.2556), 10 μg/mL loperamide ([M+H] ⁺ 477.2309), 10 μg/mL leucine ekephalin ([M+H] ⁺ 556.2771), 100 μg/mL amiodarone ([M+H] ⁺ 646.0315), 10 μg/mL
Test Compounds	28 compounds, mass range 260 – 652 amu 333 μM in 1:1 CH ₃ OH:H ₂ O

Automatic Calibration:

- Spectral accuracy calibration run daily off-hours at night unattended using QC function within MassLynx™ OALogin software.
- Calibration based on mass accuracy and spectral isotope pattern of each calibrant using MassWorks™ software (>98% SA, <5mDa error).
- Pass/fail calibration report emailed to system administrator.

Sample Analysis:

- Chemist selects designated method within MassLynx™ OALogin.
- Chemist inputs sample information, including (if known) target molecular formula, list of expected elements, and/or target ion. See example login page:
- Sample analyzed automatically by MassWorks™ when run is complete.
- Ions evaluated for mass accuracy and spectral isotope pattern to determine molecular formulas.
- PDF report containing MassWorks™ results emailed to chemist, in addition to any customary MassLynx™ report.

Customizable Parameters:

- Specified by lab system administrator within configuration file.
- Calibration specifications: compounds, SA, mDa error.
- Formula confirmation specification: SA, mDa error.
- Unknown compound formula specifications: default elements, and element limits (when no elements specified by user), C_{0-x}H_{0-x}N_{0-x}O_{0-x}S₀₋₂P₀₋₁Cl₀₋₃Br₀₋₃F₀₋₃, X based on compound mass.
 - Number of sample peaks to analyze.
 - Number of formulas to propose for each ion.

Conclusions

- Automated molecular formula determination is robust and easy to use in an open-access environment using a single quadrupole instrument without the need for high resolution instrumentation (TOF, Orbitrap).
- Software is highly successful (96%) at confirming known molecular formulas.
- For unknown compounds, specification of target mass greatly increases the likelihood that the correct formula is identified.
- For unknown compounds, chemist knowledge of likely elements present helps achieve optimum results.
- Automated open access formula determination has been successfully used in discovery chemistry project support.

Acknowledgments

AstraZeneca Infection and Cancer Chemistry