

# Investigation of Self-Calibrated, Fine-Isotope-Line Fit Scoring in Combined HPLC-Fourier Transform Orbital Trapping MS for Candidate Formulae Elimination

Yongdong Wang and Ming Gu, Cerno Bioscience, Norwalk, CT

Solving Unknown's Formula – Mass Accuracy Important

BUT – in FTMS At R=250,000 – Fine A+1, A+2, ... Isotope "Fingerprint" and Fine Structures Reveal Atom Content too


SO - Comparison of Theoretical and Observed Data for Fit Scores when Sorting Atomic Composition "Hit List"

Here: Application of the sCLIPS processing to Data From **Combined HPLC-FT-Orbital Trapping MS** is Shown

AND – The Problem of Changing Fine Isotope Structures during HPLC Elution Is Examined

EXAMPLE OF [A+2] ISOTOPE STRUCTURE OF "UNKNOWN"

WHICH ANSWER IS BETTER CURVE FIT? CALCULATE IT AND RANK!  
VERY RAPID EVEN FOR LONG LISTS OF "HITS"



- **Determine** the actual line shape of the monoisotope line, that is, the true instrument function

- **Transform** the actual instrument line function into a mathematically defined function, e.g. , a perfect Gaussian


- **Reapply** this transformation to the higher A+1, A+2, ... isotope lines (fine structures)

- **Compute** theoretical mass spec for each candidate on the elemental composition "hit" list derived from mass accuracy, double bond equivalents, and other parameters



- **Compare** each theoretical mass spec to the peak shape calibrated actual MS data – Fitting Residual (RMSE) & Spectral Accuracy

- **Combine** the information on mass error (magnitude and sign) with RMSE or Spectral Accuracy to arrive at correct formula

→ **self-Calibrated Line shape Isotope Profile Search (sCLIPS)**




Analyzing scan-by-scan mass spectrum with sCLIPS across the chromatographic peak



Reveals strong scan-to-scan dependency of both spectral and mass accuracy, with **TWO** spectral accuracy "sweet spots" located near but above the bottoms of the chromatographic peak.

Sweet Spot (RT=12.67Min)  
Spectral Accuracy: 97.4%

Apex (RT=12.60Min)  
Spectral Accuracy: 90.0%



**36 hits.** Correct formula as

- 21<sup>st</sup> mass accuracy hit
- 1<sup>st</sup> spectral accuracy hit

**42 hits.** Correct formula as

- 9<sup>th</sup> mass accuracy hit
- 3<sup>rd</sup> spectral accuracy hit

**Wide Open but Reasonable sCLIPS Search Elements:**  
C [1-25], H [0-46], N [0-13], O [0-13], S [0-2], F [0-11], P [0-6], I [0-2]  
Other generous conditions (4ppm, DBE etc)

Use of sCLIPS in HPLC / **Orbitrap FTMS**

ID and Ranking of Operational Qualification (OQ) Mixture

Formula – M/Z (nom)	Rank	Spectral Accuracy	Rejected List	*Reasons-Errors	
C <sub>7</sub> H <sub>11</sub> N <sub>4</sub> O <sub>2</sub> S	215	#1	92.70	#2-#3	Error sign / Magnitude
C <sub>28</sub> H <sub>38</sub> N <sub>8</sub> O <sub>7</sub>	556	#1	95.86	#2-7	Error sign / Mag. / [A+2]
C <sub>13</sub> H <sub>11</sub> N <sub>3</sub> O <sub>3</sub>	381	#1	90.90	#2-4	Error sign / Magnitude
C <sub>37</sub> H <sub>68</sub> NO <sub>13</sub>	741	#1	95.30	#2-4	Error sign / Magnitude
C <sub>12</sub> H <sub>15</sub> N <sub>4</sub> O <sub>2</sub> S	311	#1	92.70	#2-4	Error sign / Mag. / [A+2]
C <sub>23</sub> H <sub>21</sub> N <sub>2</sub> O <sub>9</sub>	609	#1	93.60	#2-7	Error sign / Mag. / [A+2]
C <sub>22</sub> H <sub>22</sub> NO <sub>2</sub>	472	#1	93.53	#2-5	Error sign / Mag. / [A+2]

*Orbitrap Errors are SYSTEMATIC and PRECISE.* See Strife, R.J. 2011 "Orbitrap High Resolution Applications," Ch. 4, in Characterization of Impurities and Degradants Using Mass Spectrometry, B. Pramanik, M.S. Lee and G. Chen (ed.), John Wiley & sons, Hoboken, NJ 2011  
Mathematically possible atom limits, restricted atom set (orthogonal information), Even electron restriction, RDB and ppm windows "wide"

**Synthetic Impurities Industrial Application**

HPLC / **Orbitrap FTMS**

sCLIPS Used in the Identification of 8 Synthetic Impurities

Formulas contain C, H, N, O, S, Na – C<sub>30</sub>-C<sub>44</sub>, range – MW 500-1000 Da

7/8 Compounds Ranked #1- #3 With Others Rejected By Error Sign / Magnitude

1/8 Ranked #4 - Others Rejected By Error Sign / Magnitude and [M+2] Inspection

## Findings...

- Measurement **variations** on all isotopes, A+1, A+2, ... All under AGC!
- "Sweet spots" do exist and very useful!
- Sign/magnitude of mass errors useful as well.
- Spectral Accuracy ≥ 98.0%  
→ **Unique** Formula!
- More trapping experiment {AGC thresholds, m/z}?

## Acknowledgement

1. Dr. Robert Strife from P&G for the experimental data and many fruitful discussions.
2. Dr. Christine Gu from Genentech for interesting discussions.

## References

1. Wang, Y.; Gu, M. Anal. Chem. **2010**, *82*, 7055.
2. Blake, S. L. et al, J. Am. Soc. Mass Spectrom., **2011**, *22*, 2269.
3. Gorshkov, M. V. et al, Rapid Comm. Mass. Spectrom., **2012**, *15*, 1711.
4. Kuehl, D., Am. Lab. Online, January **2008**.